



## T-100 MRI Conditional Patient Trolley



MRI CONDITIONAL  
TO 3 T.0 TESLA

### Troyka Med Inc

IOSB Sitesi 2284. Cad No: 48 ANKARA-TURKIYE

Tel: +90 5325128580- info@troykamed.com – [www.troykamed.com](http://www.troykamed.com)


T-100 MR Stretcher Op. Man. Pub. 022018 Pub. No: 0 Rev:0



ISO 13485:2016







## General Specifications

- The T-100 MR Trolley is designed to be MR Conditional and is manufactured and tested to the highest standards.
- The T-100 is constructed from 2 main materials 304 stainless steel and aluminum.
- This heavy-duty, all stainless-steel frames, non-magnetic stretcher is an ideal unit for MRI patient transport. Finished in attractive white powder coat, it comes complete with matters, restraining strips and retractable guard rails.
- Although light in weight for easy, controlled mobility, the MRI Stretcher accommodates up to 160 Kg patient safely and comfortably.
- It rolls on 5" ball-bearing swivel casters that lock firmly to secure the stretcher in place when positioned.



## Technical Specifications

<b>Manufacturer</b>	Troyka Med Inc
<b>Model:</b>	T-100
<b>Dimensions</b>	60x190x80
<b>Color</b>	9010 Pure white
<b>Main construction material:</b>	Stainless steel
<b>Guard rails</b>	Plastic Injection and aluminum
<b>Wight capacity</b>	160 Kg $\pm$ 10
<b>Weight</b>	48 kg
<b>Origin</b>	Turkey

	Protect from direct sunlight		Protect from rain and humidity		MRI conditional up to 3 Tesla
	Read user manual		Manufacturer		ISO13485 ISO 9001 Class I medical devises

