

## MR Uyumlu Pulse Oksimetre (SP02)

### MRI Compatible Pulse Oximeter (SP02)

- Monitör ve parmak probu arasındaki veri aktarımı kablosuz olarak sağlanmaktadır.
- Monitör MR odası içerisinde istenilen yere yerleştirilebilir.
- Monitör ve parmak probu MR uyumlu bataryaya sahiptir.
- Sistemin kontrolü, dokunmatik ekranla sağlanmaktadır.
- Alt ve üst alarm limitleri ayarlanabilmekte, sesli ve görsel olarak uyarı vermektedir.
- Yetişkin, pediyatrik ve yeni doğan için özel problemleri mevcuttur.
- Herhangi bir aksesuar veya sarf malzemesine ihtiyaç duymamaktadır.
- 3 Tesla'ya kadar MR uyumluluk sertifikasına sahiptir.
- Wireless data transfer between monitor and finger probe.
- Monitor can be placed anywhere in MRI room.
- Integrated MRI compatible battery in monitor and finger probe.
- System is controlled on touch screen.
- Adjustable light and voice alarms.
- Adult, pediatric, neonatal probe options.
- No need for extra accessories.
- MRI compatibility certificate up to 3 Tesla.





...sing Plug  
0997 mm

SPEED

START/STOP

TROYKA  
TROYKAMED SCREEN 1  
SPO2 %  
96  
PULSE  
87