



MRI LED Monitor



MRI CONDITIONAL
TO 3 T.0 TESLA

Troyka Med Inc
IOSB Sitesi 2284. Cad No: 48 ANKARA-TURKIYE
Tel: +90 312 265 00 96- info@troykamed.com – www.troykamed.com
T-32 MRI LED Monitor



ISO 13485:2016

MRI LED Monitor

Multipurpose Usage

This brand new, stunning 32" HD MRI LED monitor was designed to provide an optimal MRI compatible monitor that satisfies the needs of both clinical and advanced scientific applications.

Optimal Design


With its slim design, high definition display and superior image quality, the MRI LED Monitor is an optimal choice for an easy to use alternative to conventional projectors or goggle-based image delivery systems in MRI room.

Flexible Positioning

The low weight and height adjustable mobile foot stand allow easy positioning of the monitor anywhere in the MRI room. The monitor can also be wall mounted.

Patient Entertainment System

A wireless tablet with the real-time media integration playback offers patient amazing experience during the scan. The entertainment software includes over 250 animation, videos and images that patient can choose to watch during the MRI scan. Once the patient chooses the video, the video is transmitted through the wireless communication network to the monitor and the audio signal goes to the MRI sound system for the complete patient experience.



Troyka Med Inc

IOSB Sitesi 2284. Cad No: 48 ANKARA-TURKIYE

Tel: +90 312 265 00 96- info@troykamed.com – www.troykamed.com

T-32 MRI LED Monitor



ISO 13485:2016