



EMERGENCY

SIRIO S2/T

Emergency portable ventilator



Easy ventilation, advanced technology

At the 4th edition of a project started on 1979, the new **Sirio S2/T** combines excellent flexibility, friendly user and advanced technology.

The newest of Siare has been conceived expressly for patient's breathing assistance in emergency health care during transport.

Compact, light and equipped with long lasting internal battery granting the maximum flexibility of use.

Versatility:

Sirio S2/T offers a wide range of therapeutic options, it grants a controlled volume ventilation limited by pressure in invasive and non-invasive modality that meets the breathing requirements of adult, paediatric and neonatal patients (weight > 3 kg).

Simplicity:

The friendly user interface, cleverly designed to facilitate emergency health operators, makes of this ventilator the ideal choice for application on ambulance, helicopter and airplanes.

Performances:

Sirio S2/T grants excellent performances taking advantage of the technological benefits offered by Siare. It is equipped with several ventilation modalities: CPAP - PSV (Apnoea Back-up) - IPPV+AST - MANUAL – PEEP – Pressure Trigger – FiO₂ adjustable at 50% or 100% Built-in spirometry. Long lasting internal battery and a series of optional accessories designed for an easy mobility of the device.

MRI Compatible Version

The new **Sirio S2/T** is available also in the non-magnetic version (code 09636/RMN) for application in Magnetic Resonance environments from 1,5 to 3.0 Tesla. It is the ideal choice for intra-hospital transportation to transfer the patients from resuscitation to magnetic resonance in continuous ventilation without disconnection, ensuring greater comfort and safety to the patient.

SIRIO S2/T



Sirio S2/T MRI compatible can be used in helicopter.



Sirio S2/T MRI compatible can be used in Magnetic Resonance from 1,5 to 3,0 Tesla.

Friendly user settings

Acoustic and visual alarms

Built-in spirometer



Patient circuit with expiratory valve and proximal flow sensor.

Adjustable trigger

Long lasting internal battery (3 hours operation)




Version for transport with professional bag for ventilator and 2 Lt cylinder housing, complete with support for rapid connection/disconnection for installation in ambulance/helicopter.



Version for transport with professional stainless-steel tubular support for ventilator and 2 Lt cylinder housing, with possibility to attach a rapid disconnection plate for installation in ambulance/helicopter.

DATA FOR ORDERING



Code 09636 – SIRIO S2/T



Code 09636/RMN - SIRIO S2/T for MRI

OPTIONAL ACCESSORIES

 <p>Code A26.049048 Stand for ventilator with 5 wheels</p>	 <p>Code G00203001 Plate for rapid connection of tubular support 10G (code G00208003) on wall.</p>	 <p>Code 104631/SLU Silicone adult/paediatric patient circuit with expiratory valve and flow sensor Code 104631/PLU Disposable PVC single tube adult/paediatric patient circuits with exp. valve and flow sensor</p>
 <p>Code 05048/RMN Castor LT Non-magnetic trolley with 4 wheels and 1 tray</p>	 <p>Code G00213000 Plate for rapid connection of ventilator on ambulance flat surfaces</p>	 <p>Code A17.07000001 Proximal flow sensor</p>
 <p>Code 06750/SB Aluminium patient circuit supporting arm</p>	 <p>Code G90003012 Padded and refracting transport bag for ventilator and 2lt O2 cylinder complete with ambulance connection.</p>	 <p>Code 03000 VAPOUR 2 electronic humidifier</p>
 <p>Code 10006/2 Pin Index O2 cylinder 2 lt capacity, empty Code 10026/2 Pin Index O2 cylinder 2 lt capacity, empty for MRI Code 109006/CM Pin Index O2 pressure reducer</p>	 <p>Code G00208003 Stainless steel support with housing for ventilator and 2 lt. O2 cylinder complete with safety band and bag for accessories.</p>	 <p>Code G80063001 Long life external additional battery pack (for 12 hours additional operation)</p>

For additional accessories please refer to the price list in force

Distributed by:



siare
ENGINEERING
INTERNATIONAL GROUP s.r.l.

Via G. Pastore, 18 - Località Crespellano
40053 Valsamoggia (BO) Italy
Phone: +39 051 969802 (r.a.) - Fax: +39 051 969366
General e-mail: mail@siare.it - Internet: www.siare.it

