



# TH-120 Forehead Thermometer





## Screen backlight fever reminders



32°C~37.3°C

green  
normal



37.4°C~38°C


orange  
note!



38.1°C~39.1°C

red  
fever!





1. Surface mode Symbol
2. Ambient mode Symbol
3. Body mode Symbol
4. Time mode Symbol
5. Digital readout
6. Expression Symbol
7. Battery Symbol
8. Buzzer symbol
9. The order number
10. memory
11. Temperature °C (Celsius) / °F (Fahrenheit) Scale





"+"button

"-"button

Settings button

Open/Temperature measuring button

Streamlined handle

Battery cover





### Product parameters

- Product Name: Non contact IR Thermometer
- Product Model: LZ600
- Product Size: 152\*116\*43mm
- Product weight: 118g



